



Standard Disintegrator

Model 500-1



Disintegration of a pulp sample is defined as a mechanical treatment in water so that interlaced fibres, which were free in the pulp stock, are again separated from one another without appreciably changing their structural properties.

The Labtech Standard Disintegrator includes pivoted drive unit with a stainless steel propeller and shaft. Interlock protection permits safe operation only after the container is in position and the top section is locked in place. The lightweight polyurethane container supplied with the disintegrator will not dent or crack if dropped.

ABOUT LABTECH

Labtech Instruments inc. is a manufacturer of instruments used in the evaluation of Pulp and Paper properties in order to assess product characteristics and final quality.

CONTACT US

Labtech instruments inc.
Address : 1011, Autoroute 13, Laval (Québec), H7W 4V3 Canada
Phone number : 1+ (450) 681-9286
E-mail : info@labtech.ca
Fax : 450-681-2568



Features

- Integrated sensors for operator safety
 - Container Proximity Sensor
 - Cover Lock Sensor
- Rubberized container centering guides
- Lightweight polyurethane or aluminum plated Nickel container
- Digital counter
- Push button operation
- Supplied with propeller checking gauge
- Stainless steel propeller and shaft
- Splash cover
- Auto stop
- Complies with TAPPI T-205 & T-218, PAPTAC C.4 & C.5, ISO 5263



Ordering Info

Model: 500-1

Elect.: 110V/60Hz or 220V/50Hz

Instrument size: 12" x 21" x 17"

Weight: 60 lbs

Shipping dimensions: 27" x 17" x 25"

Shipping weight: 85 lbs

ABOUT LABTECH

Labtech Instruments inc. is a manufacturer of instruments used in the evaluation of Pulp and Paper properties in order to assess product characteristics and final quality.

CONTACT US

Labtech instruments inc.
Address : 1011, Autoroute 13, Laval (Québec), H7W 4V3 Canada
Phone number : 1+ (450) 681-9286
E-mail : info@labtech.ca
Fax : 450-681-2568



LABTECH INSTRUMENTS INC.

Consummables



ABOUT LABTECH

Labtech Instruments inc. is a manufacturer of instruments used in the evaluation of Pulp and Paper properties in order to assess product characteristics and final quality.

CONTACT US

Labtech instruments inc.
Address : 1011, Autoroute 13, Laval (Québec), H7W 4V3 Canada
Phone number : 1+ (450) 681-9286
E-mail : info@labtech.ca
Fax : 450-681-2568